

PRODUCT DATASHEET

TeleCare HAPA

Easy alarm release!

Explanation:

The TeleCare Handicap Adaptor (HAPA) is an auxiliary equipment to the portable trigger ATOM, which makes it very easy to activate the transmitter. The unit is intended to be used by elderly or handicapped persons that have difficulties to handle the ordinary button.

To activate HAPA you may use your palm, elbow, face or any other part of your body. The activation area is 70 x 50 mm and it requires very low pressure for activation. The unit may be worn as a necklace but may also be mounted on the wall.

A standard ATOM is easily placed into the HAPA and it can be removed if required. Note that the HAPA is delivered without the portable trigger ATOM!

Vantages:

- Unit size 102x50 mm
 - Activation area 70x50 mm
 - Worn as a necklace or wall mounted
 - Very easy activation
 - Used together with standard ATOM
-